

Cumulative Release Notes for DSQM.

Version 01.00

Features.

This updates DSQ to DSQM as follows:

1. Supports more than one hospital
2. Minor enhancements to HC21 forms

Installation.

1. Copy the three files to appropriate folder. Refer to the user guide.
2. Under Advance>Relink Databases relink the data and episodes files to the front end.
3. Start Blue Chip and point it to the DSQM.mdb file from Setup>Practice>Links>Day Surgery.

Modified File List

File Name	Version	Modified from last release
DSQ.mdb	01.01	No
DSQ Data.mde	01.01	No
DSQ Episodes.mde	01.01	No

Version 02.00

Features.

This updates DSQM as follows:

1. ICD10-AM tables updated to version 12
2. HCP based extracts from 1st July 2022 updated to report correct coding version

Installation.

1. Make backups of **EVERYTHING**
2. Replace DSQM.mdb and DSQM Data.mde.
3. Start DSQM and Under Advanced relink the data and episode databases

Modified File List

File Name	Version	Modified from last release
DSQM.mdb	02.01	Yes
DSQM Data.mde	02.01	Yes
DSQM Episodes.mde	01.02	Yes. New installs only